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Stand: 2021-04

Product description

Orthodontic resin, autopolymer. Acrylic for the fabrication of orthodontic devices using the dough technique (**orthoXin®**) or the “Salt and Pepper” technique (**orthoXin® ST**).

Composition

2-component resin system consisting of powder and liquid. Main powder components: methyl methacrylate copolymers. Main liquid components: methyl methacrylate, dimethacrylate.

Indications

- Removable orthodontic devices, short-term use
- Removable orthodontic devices, cumulated long-term use
- Fixed orthodontic devices
- Occlusal splints
- Bleaching, snoring and grinding splints
- Drilling templates for implantology
- Additions and repairs

Contraindications

Patient contact with the unpolymerized acrylic or its individual components is contraindicated. Do not use in case of known allergy to one or more components of the product.

Model fabrication

A vacuum mixer and base system are required for fabricating working models. The plaster to be used should be approved by the manufacturer for use in the orthodontic technique, or acrylic technique. The use of Class III dental plaster is recommended.

Preparation

It is essential to soak the model until air bubbles are no longer released. This is particularly important to ensure that no air is pressed into the acrylic from the plaster and produces porosity.

Separating

Apply a thin layer of alginate-based separating solution 1-2 times to the model. The **AcrylX® Xeparator** always shows reliable results.

Important: If separating twice, the first layer of separating agent should not be completely dry before applying the second layer. This guarantees an optimum bond between the layers of separating agent. The separating agent must, however, be completely dry before contact with acrylic to prevent any white discoloration on the surface of the acrylic. It is essential to follow the instructions for use of the separating solution!

Processing in the dough technique - orthoXin®

Mixing ratio 10:5

The recommended mixing ratio is 10g powder and 5g liquid. Place the liquid in a mixing beaker, add the corresponding amount of powder as quickly as possible and mix thoroughly to a homogeneous consistency using a spatula. Ensure there are no bubbles. Free dosing is possible, but a guarantee of the physical material properties is only given if the components are weighed out. With free dosing, add as much powder to the liquid until the liquid is completely absorbed.

With the beginning of the plastic phase the dough may be applied on to the plaster model.

Processing in the „Salt and Pepper“ technique - orthoXin® ST

The polymer and monomer are applied alternately to the plaster model, with the polymer first being sprinkled on and then wetted with the monomer. This process must be repeated until the desired layer thickness has been achieved. It is particularly important to note that only as much liquid may be applied as can be absorbed by the powder. Too much monomer leads to deliquescence of the material. Before polymerizing the device, the surface must be wetted again with the liquid and smoothed.

Note

When using very opaque colors, particularly the “black and white” version, ensure that expansion screws and wire retainers are not exposed when trimming and cannot be damaged. Motif inserts and shade modifiers can also lose their effectiveness. The polymer can generally be thinned using transparent material.

Polymerization

The acrylic should be polymerized under a pressure of 2-3 bar and a water temperature of 40°C [104°F] for approx. 10 minutes.

After completion

Until incorporation the device should be stored in water. It is always recommended to water the device for 24 hours to achieve further reduction of the residual monomer content (orthoXin® 3,4% / orthoXin® ST 4,4% according to ISO 20795-2) and thus prevent any taste irritation.

Storage and Safety Measures

Do not use the material after the expiration date. Store at temperatures not exceeding 25°C [77°F]. Protect from sunlight. Close containers immediately after use. Keep away from heat/sparks/open flames/hot surfaces. No smoking. Dispose of contents/ container in accordance with local disposal regulations.

Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation. Avoid breathing dust/fume/gas/mist/vapours/spray. Wear protective gloves/protective clothing/eye protection/ face protection. IF ON SKIN: Wash with plenty of soap and water.

Feedback

Please direct any feedback regarding this product to **info@acrylx.com**. All serious incidents that have occurred in connection with the product must be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is located.

Colors

The material is available in colorless, but can be customized with ready-colored monomers or color concentrates.

Delivery forms

Powder: 1kg, 5kg, 10kg
Liquid: 500ml, 5L
Additional sizes upon request

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® = registered trademark of AcrylX® GmbH